


Set-up Instructions

- 1 Fill out the PATIENT DEMOGRAPHICS section on the requisition form and have the patient sign the PATIENT CONSENT at the bottom, circled in the example.
- 2 Using a pen, write the name and date of birth of the patient on the specimen tube.

If ordering a test electronically through Copia, a requisition will print off for the patient to sign.



PATIENT DEMOGRAPHICS - Please print clearly to ensure proper delivery of your results

First Name: _____ Last Name: _____

Gender: Male Female Other Date of Birth: _____ Phone: _____

Street Address: _____ City: _____ State: _____ Zip: _____

SELECT TEST

PCR COVID-19 TEST OPTIONS (Price for test only available)

<input type="checkbox"/> PCR 24-Hour Results	\$150.00
<input type="checkbox"/> PCR 2-Hour Results	\$200.00



PATIENT CONSENT

I have read the above information and understand the risks and benefits of this test. I hereby consent to the collection, use, and disclosure of my personal information for the purposes of this test. I understand that my information may be shared with other healthcare providers for the purpose of my care. I understand that my information may be used for research and quality improvement purposes. I understand that my information may be disclosed to other parties for the purpose of my care. I understand that my information may be disclosed to other parties for the purpose of my care. I understand that my information may be disclosed to other parties for the purpose of my care.


Signature: _____ Date: _____

Collection Instructions

- 3 Remove the swab from the package. Do not touch the soft end with your hands or anything else.
- 5 Slowly rotate the swab, gently pressing against the inside of your nostril for a total of 10 times. Get as much nasal discharge as possible on the soft end of the swab. Remove the swab.




- 4 Insert the entire soft end of the swab into your nostril about 1 inch.
- 6 Using the same swab, repeat steps 4 & 5 in your other nostril. When complete remove swab the from nose and don't allow it to touch your hand or any surfaces.



Prepare Specimen

- 7 Holding the sterile tube upright, being careful not to spill any of the solution, open the tube. Place the swab in the tube and snap off the end of the swab at the break line, so that it fits in the tube. Place the cap on the tube and screw down tightly to prevent leakage.
- 9 Place the filled out and signed requisition form in the biohazard bag and ship to our lab.



- 8 Place the tube containing the swab in the provided biohazard bag.



Results

- 10 Results are available in our patient portal within 24-hours of sample receipt. Patients will be notified via email and text message when results are ready.